

Thermo Scientific™ Model 5030i SHARP Synchronized Hybrid Ambient Real-time Particulate Monitor

Spare Parts Recommendations

Part Number	Description	Annual Maintenance	Age of Instrument or Fleet (Years)					
			1	2	3	4	5	
8630	Filter Guard Assembly (w/ foam)	X	X	X	X	X	X	X
109612-00	Fuse, 250 VAC, 10.0 Amp, SlowBlow (for 110 VAC and 115 VAC models)	X	X	X	X	X	X	X
109613-00	Fuse, 250 VAC, 5.0 Amp, SlowBlow (for 220 VAC and 240 VAC models)	X	X	X	X	X	X	X
59-008630	Pump Rebuild Kit	X	X	X	X	X	X	X
106443-00	Filter Tape	X	X	X	X	X	X	X
24-000483	Pump Tubing, External 3/8-inch	X		X	X	X	X	X
106994-00	Tubing, Internal 3/8 inch			X				
107000-00	Fitting, 5/8 union w/Nylon Ferrules			X				
FH125C14	Calibration Kit (Cal Foils)			X				
10-001403	External Pump Assembly				X			
101688-00	Ambient Temperature Thermistor (with connector)					X		
106946-00	Pressure Board Assembly (w/tubing assembly)					X		
100907-00	Fan, 24 VDC							X
101055-00	AC Receptacle Assembly							X
101681-00	Power Supply Assembly, 24 VDC (with base plate and screws)							X
106470-00	Proportional Valve Assembly							X
106535-00	Heater Tube Assembly, 3 feet							X
106536-00	RH Temperature Assembly							X
106923-00	Pressure Valves							X
109976-01	5014i/5030i Enclosure with Heater & A/C - 120V 60hz							X
109976-02	5014i/5030i Enclosure with Heater & A/C - 220/240V 50hz							X
110071-01	5014i/5030i Enclosure with Heater only - 120V 60hz							X
110071-02	5014i/5030i Enclosure with Heater only - 220/240V 50hz							X

© 2014 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details. This product is manufactured in a plant whose quality management system is ISO 9001 certified.

USA
27 Forge Parkway
Franklin, MA 02038
Ph: (866) 282-0430
Fax: (508) 520-1460
customerservice.aq@thermofisher.com

India
C/327, TTC Industrial Area
MIDC Pawane
New Mumbai 400 705, India
Ph: +91 22 4157 8800
india@thermofisher.com

China
+Units 702-715, 7th Floor
Tower West, Yonghe
Beijing, China 100007
+86 10 84193588
info.eid.china@thermofisher.com

Europe
Takkebijsters 1
Breda Netherlands 4801EB
+31 765795641
info.aq.breda@thermofisher.com